




Fig. 1




Hexammine Cobalt




Chloropentammine Cobalt




Spermidine



Manganese Chloride



Spermine



(D or L) Polylysine

Figure 1. Structures of common compaction agents

Fig. 2




Figure 2. Precipitation by spermidine of 40 µg/ml pBGS19Luxwt or Baker's yeast RNA in 10 mM Tris buffer at pH 8.0 with and without 600 mM NaCl. Error bars are +/- one standard deviation.

Fig. 3




Figure 3. 1% agarose gel that traces the large-scale purification of pBGS19luxwt. Lane 1 is a supercoiled plasmid ladder from Gibco; Lane 2 is the preparation after Celite filtration, isopropanol precipitation, and resuspension; Lane 3 is the supernatant after LiCl precipitation; Lane 4 is the supernatant of the compaction precipitation; Lane 5 is the resuspended pellet of the compaction precipitation; Lane 6 is a 10X loading of the material in Lane 5; Lane 7 is after a Q sepharose anion exchange column (Fig. 5, bottom, Peak 5); Lane 8 is a 10X loading of Lane 7 and Lane 9 is pBGS19Luxwt plasmid DNA separated using the mini-prep procedure.

Fig. 4




Figure 4. FPLC anion-exchange separation of pBGS19Luxwt of an alkaline lysate after isopropanol and LiCl precipitation. Top: NaCl gradient; Middle: with no previous compaction precipitation step; Bottom: identical separation after a previous compaction precipitation step.

Fig. 5




Figure 5. Summary of a selective precipitation-based noncolumn DNA purification as detailed in example 1.

Fig. 6




Figure 6. 3% biogel (from Bio101 Inc.) of *V. proteolyticus* RNA purified by Example 9. Lane 1 is the Ambion RNA Century Plus Size Markers; Lane 2 is the lysate after BPER addition, spermidine addition and centrifugation; Lane 3 is the supernatant of the 4 mM hexamine cobalt precipitation; and Lane 4 is the RNA pelleted in the hexamine cobalt precipitation but before any column separation.

Fig. 7




Figure 7. FPLC chromatogram of *V. proteolyticus* RNA on and 25 ml high performance Q Sepharose anion-exchange column (Pharmacia). The gradient was ran over 12 column volumes from 0.30 M NaCl to 0.57 M NaCl in a column buffer of 20 mM bis-tris propane and 20 mM EDTA at pH 6.9.

Fig. 8




Figure 8. FPLC chromatogram of pCP3X3 aRNA containing *E. coli* strain JM109 on and 25 ml high performance Q Sepharose anion-exchange column (Pharmacia). The gradient was ran over 12 column volumes from 0.37 M NaCl to 0.57 M NaCl in a column buffer of 20 mM bis-tris propane and 20 mM EDTA at pH 6.9.

Fig. 9





Figure 9. FPLC chromatogram of selective precipitation purified β ribozyme on and 25 ml high performance Q Sepharose anion-exchange column (Pharmacia). The gradient was ran over 12 column volumes from 0.37 M NaCl to 0.7 M NaCl in a column buffer of 10 mM bis-tris propane and 2 mM EDTA at pH 6.9.

Drawing
10/1

Description of Drawing 10B

10/1
24/10/01



KNA isolation kit based on a hexamine cobalt precipitation process

- Lysis/DNA precipitation solution (e.g. 50% HFEP with 2.5 mM Spermidine or 1% Rnaji 58 with 2.5 mM spermidine in 10 mM bis tris propane at pH 7). This may also be extended for use with plant cells and other eukaryotic cells with the possibilities of homogenisation, other lysis solutions, and breaking the lysis from the spermidine DNA removal. Thus, there would be a lysis solution and a separate DNA precipitation solution.
- Hexamine cobalt precipitation solution (e.g. 7 mM hexamine cobalt for a total KNA precipitation or 4 mM hexamine cobalt in 10 mM bis tris propane at pH 6.9 for a high molecular weight precipitation).
- Optional: a second hexamine cobalt precipitation solution to bring down low molecular weight RNA not precipitated when solution B was used (e.g. 20 mM hexamine cobalt in 10 mM bis tris propane at pH 6.9).
- Stripping solution (e.g. 50% isopropyl alcohol with 3 M Urea, 300 mM NaCl, 25 mM EDTA).
- 70% EtOH wash^{1/2}
- Final Resuspension solution (e.g. High purity TE which is 10 mM Tris and 1 mM EDTA at pH 8.0).^{1/2}

^{1/2}User may have to add enough alcohol to bring the solution to the proper percentage of alcohol.

^{1/2}Optional kit components that could be provided by end user.

Drawing 10/1
10/1

Description of Drawing 10A

Plasmid DNA separation kit based on competition agent precipitation technology. In the box a total of 8 solutions should be included two of which are optional. Solutions will include the 3 common alkaline lysis solutions, a low ionic strength resuspension buffer, a competition agent precipitation solution, and a stripping solution, and optionally a 70% ethanol wash solution and a final resuspension solution.


Description of each bottle are as follows:

- Alkaline lysis solution I (e.g. 25 mM Tris and 10 mM EDTA at pH 8.0)
- Alkaline lysis solution II (e.g. 1% Sodium Dodecyl Sulfate (SDS) 0.2 N NaOH)
- Alkaline lysis solution III (e.g. 3 M KAc at pH 5.5)
- Resuspension solution (e.g. 10 mM Tris at pH 8.0)
- Competition agent precipitation solution (e.g. 2 mM Spermidine 3HCl and 10 mM Tris at pH 8.0)
- Competition agent stripping solution (e.g. 50% EtOH, 300 mM NaCl, 12.5 mM EDTA)
- 70% EtOH wash^{1/2}
- Final Resuspension solution (e.g. High purity TE which is 10 mM Tris and 1 mM EDTA at pH 8.0).^{1/2}

^{1/2}User may have to add enough EtOH to bring the solution to the proper percentage of EtOH.

^{1/2}Optional kit components that could be provided by end user.

Also, solutions D and E can be combined to form a resuspension/competition agent precipitation solution.






Fig. 11




Fig. 12






Fig. 13



5


FIG. 15



009MUS20010418

FIG.16

5




10

15

FIG.17

1 2



5

FIG.18

